


PHYSIO FLOW™

Noninvasive Cardiac Output

PHYSIO FLOW **ENDURO**

ULTRALIGHT AND ACCURATE



NEW : PhysioFlow Parameters

In a holter size device

- Stroke Volume
- Cardiac Output / Index
- Contractility
- Ventricular Ejection Time
- Ejection Fraction (estimate)
- Early Diastolic Function ratio
- Systemic Vascular Resistance
- Left Cardiac Work index

For Multiple applications

- Cardiology
- Stress testing
- Physiology and sports medicine
- Military and aerospace medicine

Holter Monitor Features

- Small Size 115 x 85 x 18 mm
- Light Weight : Less than 200g (with batteries)
- Standard AA batteries or rechargeable, 6 hours autonomy
- 24 hours MMC memory
- USB or BlueTooth wireless memory download
- Real time Wireless monitoring using BlueTooth
- Works with PhysioFlow PF106 MS Windows™ XP based Software for display, data analysis, and storage

Features of device may change. Governmental approvals pending

